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U-LAB

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Instructions for use

U-LAB

Chemistry test card for the rapid determination of Nitrite and Leukocytes in urine. The test is an in vitro diagnostic test intended to be used by lay people for self-testing of urine.

SUMMARY AND EXPLANATION OF THE TEST

U-Lab is used as a screening test for the early identification of urinary tract infections for persons at risk. The test is non invasive and simple to use. The card shows a test chamber containing dry chemistry reagent strips for Nitrite, Leukocytes and sample detection. In the event of a possible pathological finding (= positive test), a color change occurs in the respective test area. The color reaction shows a "positive or negative result" for Nitrite and Leukocytes in urine.

PRINCIPLES OF THE TEST

Nitrite: The test indirectly detects the presence of bacteria that reduce Nitrate, usually present in urine, into Nitrite. The test pad for Nitrite is impregnated with an amine and a coupling reagent. Nitrite present in the urine diazotize the amine. The subsequent coupling reaction produces a pink color.


Leukocytes: The Leukocytes excreted in the urine are almost exclusively granulocytes, whose esterase activity is detected in the correspondent test chamber. The test zone contains an indoxyl ester and a diazonium salt. Granulocyte esterases split the ester, and as a result, the free indoxyl can react with the diazonium salt to produce a violet color.

TEST PROCEDURE

Each test kit provides 2 or 6 tests. The U-Lab test kit consists of 2 or 6 individually packed test cards, 2 or 6 plastic containers and the instructions for use leaflet. A pictogram on the back of each box explains how to perform the test. However, this pictogram does not replace the instructions for use which is the only legal valid document explaining its usage.

Method:

1. Remove one plastic container from the package.
2. Collect your first morning urine into the plastic container.
3. Take the test card out of its individual package, drop the card into the urine and leave it in the container for 5 minutes.
4. After 5 minutes, take the card out of the urine and compare the color reaction in the test chamber with the reference color chart.
5. Color changes occur already after 2-3 minutes but results are only valid after 5 minutes.
6. If there is any doubt concerning the result perform a second test the next morning.
7. The color reaction for Leukocytes and Nitrite remains stable for 24 hours.

I.

II.

III.

INTERPRETATION OF RESULTS

Results are determined visually by direct comparison of the reacted test fields in the test chamber, with the reference color chart above each test field.

Nitrite indicator: pink color = test is positive


Leukocytes indicator: violet color = test is positive

Sample volume indicator: orange color = sufficient sample in the chamber


Non valid test:

not enough urine in the chamber

Valid test:

Nitrite negative, Leukocytes negative

Valid test:

Nitrite positive, Leukocytes negative

Valid test:

Nitrite negative, Leukocytes positive

Valid test:

Nitrite positive, Leukocytes positive

REAGENTS

Valid: methyl red 2,0%, bromothymol blue 10,0%
Leukocytes: carboxylic acid ester 0,4%, diazonium salt 0,2%
Nitrite: tetrahydrobenzo(h)quinolin-3ol 1,5%, sulfanic acid 1,9%

LIMITATIONS OF THE TEST

1. The test is intended for screening only. All positive results (for Nitrite and/or Leukocytes) should be confirmed by a second test.
2. If the second test is again positive (for Nitrite and/or Leukocytes), you are advised to see your medical doctor.
3. Diagnostic or therapeutic decisions should not be based on any single result or method but based on an in depth examination by your medical doctor.

Nitrite: Even in the presence of bacteria, the test can show a negative result because of bacteria that do not have the capacity to reduce Nitrate into Nitrite. False-negative responses, can be caused by the following: Antibiotic therapy, low-nitrate diets, strong diuresis, high levels of ascorbic acid, high specific gravity or insufficient urinary retention in the bladder. False-positive responses can be caused by dyes excreted in the urine (e.g. pyridium, red beets).

Leukocytes: False-positive reactions may be caused by formaldehyde (Preservation). Protein concentrations > 500 mg/dl may diminish the color response. Bacteria, trichomonads and red blood cells do not, however, react with the test. High daily doses of cephalexin or gentamicin as well as very high glucose concentrations or high specific gravity may diminish the color response.

EXPECTED VALUES

Nitrite: This test responds to urinary Nitrite levels from 0,05 to 0,1 mg/dl by developing a pink color. Normally, no Nitrite are detectable in urine and therefore, any positive response should be considered indicative of a urinary tract infection of bacterial origin. However, a negative result does not rule out an infection of the urinary tract.

Leukocytes: Normal urine often contains 0 - 2 Leukocytes per microliter urine. The practical sensitivity of this test has been set at 10 - 20 Leukocytes per microliter urine. Partial cytolysis intensifies the color response, particularly in the region of the maximum analytical sensitivity.

WARNINGS AND PRECAUTIONS

1. All unused cards must remain in their original package.
2. Do not open the original package if you are not using the test card immediately.

STORAGE AND STABILITY

Store the package at 2°C to 25°C under dry conditions. Unopened test cards capped in the original box, are stable until the expiration date.

0483 = declared conform

= for in vitro diagnostic use only

= this card is for a single use

REF = article-number

= manufacturing date

= store in a dry place, storage temperature 2 - 25°C

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